

Pfu Mix (2x)

Catalog No: P1092 (5x1ml)

Contents: 2xPfu Mix 5x1 ml
 Nuclease-free water 5x1 ml

Store at -20°C

In total 10 vials.

Description

Pfu Mix(2x) is a premixed, ready-to-use solution containing Pfu DNA Polymerase, dNTPs, MgSO₄ and Reaction Buffer at optimal concentrations for efficient amplification of DNA templates by PCR. To prepare the final PCR, only primers and template DNA are added. Pfu Mix contributes to highly reproducible PCR by reducing the risk of pipetting errors, miscalculation and contamination. It also contributes to higher sensitivity by adding intensifier and optimizer.

Pfu DNA polymerase, derived from the hyperthermophilic archae *Pyrococcus furiosus*, has been shown to exhibit superior thermostability and proofreading properties compared to other thermostable polymerase. Unlike Taq DNA polymerase, highly thermostable Pfu DNA polymerase possesses 3' to 5' exonuclease proofreading activity that enables the polymerase to correct nucleotide-misincorporation errors. This means that Pfu DNA polymerase-generated PCR fragments will have fewer errors than Taq-generated PCR inserts. Using Pfu DNA polymerase in your PCR reactions results in blunt-ended PCR products, which are ideal for cloning into blunt-ended vectors. Pfu DNA polymerase is superior for techniques that require high-fidelity DNA synthesis.

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

Web: b-kl.eu
e-mail: b-kl@b-kl.eu

PRODUCT USE LIMITATION.

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Applications

- High fidelity PCR
- Routine PCR with high reproducibility
- Generation of PCR products for TA cloning
- Site-directed mutagenesis

Feature

- Convenient –Pfu DNA Polymerase in a ready-to-use Mix.
- High yields of PCR products with minimal optimization.
- Fast -saves time due to reduced number of pipetting steps.
- Reproducible -lower contamination and pipetting error risk.

Composition of the 2x Pfu Mix

0.15U/ul Pfu DNA polymerase, 2xPfu Buffer, 0.4mM dNTPs, 4mM MgSO₄, 0.02% bromophenol blue.

Pfu mix buffer is a proprietary formulation optimized for robust performance in PCR.

Basic PCR Protocol

All solutions should be thawed on ice, gently vortexed and briefly centrifuged.

1. Add in a thin walled PCR tube on ice:

For a total 50µl reaction volume

Component of sample	Volume	Final concentration
Pfu Mix (2X)	25 µl	1X
Forward Primer	variable	0.1-1 µM
Reverse Primer	variable	0.1-1 µM
Template DNA	variable	10 pg-1 µg
Water, nuclease-free	to 50 µl	–

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

Recommendations with Template DNA in a 50µl reaction volume

Human genomic DNA	0.1 µg-1 µg
Plasmid DNA	0.5 ng-5 ng
Phage DNA	0.1 ng-10 ng
E.coli genomic DNA	10 g-100 ng

2. Gently vortex the sample and briefly centrifuge to collect all drops to the bottom of the tube.

3. Overlay the sample with mineral oil or add an appropriate amount of wax. This step may be omitted if the thermal cycler is equipped with a heated lid.

4. Perform PCR using the following thermal cycling conditions.

Initial Denaturation	94°C	3 minutes
25-35 Cycles	94°C	30 seconds
	55-68°C	30 seconds
	72°C	2-3 minutes
Final Extension	72°C	10 minutes

5. Incubate for an additional 10 min at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.

6. Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide staining. Use appropriate molecular weight standards.

Notes on cycling conditions

- Pfu Mix is for High fidelity PCR applications.
- The error rate of Pfu DNA Polymerase in PCR is 2.6×10^{-6} errors per nt per cycle; the accuracy (an inverse of error rate) an average number of correct nucleotides incorporated before making an error is 3.8×10^5 (determined according to the modified method described in)
- Pfu DNA Polymerase accepts modified nucleotides (e.g. biotin-, digoxigenin-, fluorescent-labeled nucleotides) as substrates for the DNA synthesis.
- The enzyme has no detectable reverse transcriptase activity.
- The number of PCR cycles depends on the amount of template DNA in the reaction mix and on the expected yield of the PCR product. 25-35 cycles are usually sufficient for the majority PCR reaction. Low amounts of starting

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

template may require 40 cycles.

Guidelines for preventing contamination of PCR reaction

During PCR more than 10 million copies of template DNA are generated. Therefore, care must be taken to avoid contamination with other templates and amplicons that may be present in the laboratory environment. General recommendations to lower the risk of contamination are

as follows:

- Prepare your DNA sample, set up the PCR mixture, perform thermal cycling and analyze PCR products in separate areas.
- Set up PCR mixtures in a laminar flow cabinet equipped with an UV lamp.
- Wear fresh gloves for DNA purification and reaction set up.
- Use reagent containers dedicated for PCR. Use positive displacement pipettes, or use pipette tips with aerosol filters to prepare DNA samples and perform PCR set up.
- Always perform “no template control” (NTC) reactions to check for contamination

Quality Control

The absence of endodeoxyribonucleases, exodeoxyribonucleases and ribonucleases is confirmed by appropriate quality tests. Functionally tested in amplification of a single-copy gene from human genomic DNA.

Endodeoxyribonuclease Assay

No detectable conversion of covalently closed circular DNA to a nicked DNA was observed after incubation of 25µl Pfu Mix (2X) with 1µg pBR322 DNA in 50µl for 4 hours at 37°C and 70°C.

Exodeoxyribonuclease Assay

No detectable degradation of lambda DNA-HindIII fragments was observed after incubation of 25µl Pfu Mix (2X) with 1µg digested DNA in 50µl for 4 hours at 37°C and 70°C.

Ribonuclease Assay

0% of the total radioactivity was released into trichloroacetic acid-soluble fraction after incubation of 25µl of Pfu Mix (2X) with 1µg E.coli [3H]-RNA (40000cpm/µg) in 50µl for 4 hours at 37°C and 70°C.

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.