

HSTM Kit

Catalog No.: P2061(5x1 ml)

Contents: 2X HSTM Reaction Mix 5x1 ml
 HSTM Taq DNA Polymerase (2.5U/μl) 160μl

Store at -20°C

In total 6 vials.

Description

HSTM Kit contains a HSTM Taq DNA polymerase and a 2x premixed, ready-to-use HSTM reaction Mix. The HSTM reaction mix is supplied as a 2X concentrated solution, that contains dNTPs and all other PCR components, except DNA template and primers. To prepare the final PCR reaction system, mix HSTM Taq DNA polymerase and the HSTM reaction mix, then add primer and template. The reaction mix contributes to high specificity by optimizing the system, reducing primer-dimer rate. PCR amplification sensitivity is controllable due to the amount of DNA polymerase is flexible.

HSTM Taq DNA Polymerase is a thermostable recombinant DNA polymerase derived from thermophilic bacterium *Thermus aquaticus*, its molecular weight is 94 kDa. HSTM Taq DNA Polymerase can amplify DNA target up to 5 kb. The elongation velocity is 0.9~1.2kb/min. It has 5' to 3' polymerase activity but lacks of 3' to 5' exonuclease activity, that results in a 3'-dA overhangs PCR product. All components of the HSTM Mix are at optimal concentration for efficient amplification, it contributes to highly specific incorporation of primer and template.

PRODUCT USE LIMITATION.

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

Web: b-kl.eu
e-mail: b-kl@b-kl.eu

Unit Definition

One unit is defined as the amount of the enzyme required to catalyze the incorporation of 10 nmole of dNTPs into an acid-insoluble form in 30 minutes at 70°C using hering sperm DNA as substrate.

Storage Buffer

20mM TrisCl (pH8.0), 100mM KCl, 3mM MgCl₂ 1mM DTT, 0.1% NP-40 ,0.1% Tween20, 0.2mg/ml BSA, 50% (v/v) glycerol

Composition of 2X HS™ Reaction Mix

2xHS™ PCR Buffer, 0.4mM dNTPs, 3.2mM MgSO₄, 0.02% bromophenol blue.

HS™ Reaction Mix is a proprietary formulation optimized for robust performance in PCR

Features

- **Efficiency:** simplifying the operation and saving your time
- **Reproducible:** reduce the risk of pipetting errors, miscalculation and contamination
- **Flexible:** the amount of polymerase is flexible and controllable.

Applications

- PCR amplification of DNA fragments as long as 5 kb
- DNA labeling.
- DNA sequencing.
- PCR for cloning

Basic PCR Protocol

All solutions should be thawed on ice, gently vortexed and briefly centrifuged.

1. Add in a thin walled PCR tube on ice:

For a total 50µl reaction volume

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

Component of sample	Volume	Final concentration
2X <i>HS</i> TM Reaction Mix	25 µl	1X
Forward Primer	variable	0.1-1 µM
Reverse Primer	variable	0.1-1 µM
Template DNA	variable	10 pg-1 µg
<i>HS</i> TM Taq DNA Polymerase (2.5U/µl)	0.5-1 µl	1.25-2.5U/50µl
Water, nuclease-free	to 50 µl	–

Recommendations with Template DNA in a 50µl reaction volume

Human genomic DNA	0.1 µg-1 µg
Plasmid DNA	0.5 ng-5 ng
Phage DNA	0.1 ng-10 ng
E.coli genomic DNA	10 g-100 ng

2. Gently vortex the sample and briefly centrifuge to collect all drops to the bottom of the tube. Overlay the sample with mineral oil or add an appropriate amount of wax.

This step may be omitted if the thermal cycler is equipped with a heated lid.

3. Place samples in a thermocycler and start the program
Perform 25-35 cycles of PCR amplification as follows:

Initial Denaturation	94°C	3 minutes
25-35 Cycles	94°C	30 seconds
	55-68°C	30 seconds
	72°C	1 minutes
Final Extension	72°C	10 minutes

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

4. Incubate for an additional 10 min at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.

5. Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide staining. Use appropriate molecular weight standards.

Notes on cycling conditions

- *HS*TM Taq DNA Polymerase is for High Specificity PCR applications.
- The half-life of enzyme is >40 minutes at 95°C.
- The error rate of *HS*TM Taq DNA Polymerase in PCR is 2.2×10^{-5} errors per nt per cycle; the accuracy (an inverse of the error rate) an average number of correct nucleotides incorporated before making an error, is 4.5×10^{-4} (determined according to the modified method described in).
- *HS*TM Taq DNA Polymerase accepts modified nucleotides (e.g. biotin-, digoxigenin-, fluorescent-labeled nucleotides) as substrates for the DNA synthesis.
- The number of PCR cycles depends on the amount of template DNA in the reaction mix and on the expected yield of the PCR product. 25-35 cycles are usually sufficient for the majority PCR reaction. Low amounts of starting template may require 40 cycles.

Guidelines for preventing contamination of PCR reaction

During PCR more than 10 million copies of template DNA are generated. Therefore, care must be taken to avoid contamination with other templates and amplicons that may be present in the laboratory environment. General recommendations to lower the risk of contamination are as follows:

- Prepare your DNA sample, set up the PCR mixture, perform thermal cycling and analyze PCR products in separate areas.
- Set up PCR mixtures in a laminar flow cabinet equipped with an UV lamp.
- Wear fresh gloves for DNA purification and reaction set up.
- Use reagent containers dedicated for PCR. Use positive displacement pipettes, or use pipette tips with aerosol filters to prepare DNA samples and perform PCR set up.
- Always perform “no template control” (NTC) reactions to check for contamination

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

Quality Control

The absence of endodeoxyribonucleases, exodeoxyribonucleases and ribonucleases is confirmed by appropriate quality tests. Functionally tested in amplification of a single-copy gene from human genomic DNA.

Endodeoxyribonuclease Assay

No detectable conversion of covalently closed circular DNA to a nicked DNA was observed after incubation of 10U *HS*TM Taq DNA Polymerase with 1 µg pBR322 DNA for 4 hours at 37°C and 70°C.

Exodeoxyribonuclease Assay

No detectable degradation of lambda DNA-HindIII fragments was observed after incubation of 10U *HS*TM Taq DNA Polymerase with 1 µg digested DNA for 4 hours at 37°C and 70°C.

Ribonuclease Assay

0% of the total radioactivity was released into trichloroacetic acid-soluble fraction after incubation of 10U *HS*TM Taq DNA Polymerase with 1 µg E.coli [3H]-RNA (40000cpm/µg) for 4 hours at 37°C and 70°C.

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

Web: b-kl.eu
e-mail: b-kl@b-kl.eu