
FSTM Taq DNA Polymerase

Catalog No.: P3031 (250 U)

Concentration: 5U/ μ l

Contents: FSTM Taq DNA Polymerase 50 μ l
10xFSTM PCR Buffer (Mg²⁺ Plus) 1.25ml
6xLoading Buffer 1ml

Description

FSTM Taq DNA Polymerase is the latest generation Taq-based DNA polymerase. It possesses high amplification efficiency as Taq DNA polymerase does, and fast elongation ability as KOD polymerase does, can be used in a variety of PCR. The FSTM PCR Buffer, designed for FSTM Taq DNA polymerase, can be used in fast amplification reaction. The elongation rate of FSTM Taq DNA polymerase can be as much as 3kb/min. It can shorten the amplification time by 1/3.

Unit Definition

One unit is defined as the amount of the enzyme required to catalyze the incorporation of 10 nmole of dNTPs into an acid-insoluble form in 30 minutes at 70°C using hering sperm DNA as substrate.

Storage Buffer

20mM TrisCl (pH8.0), 1 mM DTT, 0.1 mM EDTA, 100 mM KCl, 0.5% TW 20, 0.5% NP 40, 50% Glycerol

10X FSTM PCR Buffer with Mg²⁺

200 mM Tris-Cl (PH 8.8), 100 mM KCl, 16 mM (NH₄)₂SO₄, 16 mM MgSO₄, 1% Triton-X-100.

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

PRODUCT USE LIMITATION.

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Features

- Fast elongation: elongation rate can reach 3kb/min, 3 times rate of Taq DNA polymerase
- Thermostable: half-life over 40 min at 95°C incubation.
- Generates 3'-dA overhangs PCR products.

Applications

- Routine PCR
- DNA labeling
- PCR sequencing
- Generate PCR product for TA cloning

Basic PCR Protocol

The following basic protocol serves as a general guideline and a starting point for any PCR amplification. Optimal reaction conditions (incubation time and temperature, concentration of FSTM Taq DNA Polymerase, primers, Mg²⁺, and template DNA) vary and need to be optimized.

1. Add the following components to a sterile microcentrifuge tube sitting on ice:

1.1 Recommended PCR assay with FSTM PCR Buffer (Mg²⁺ plus)

| Reagent | Quantity, for 50µl reaction | Final concentration |
|---|-----------------------------|---------------------|
| Sterile deionized water | variable | - |
| 10X FS TM PCR Buffer (Mg ²⁺ plus) | 5µl | 1X |

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

| | | |
|--|------------------|----------------------|
| dNTPs (10mM each) | 1 μ l | 0.2 mM each |
| Primer I | variable | 0.4-1 μ M |
| Primer II | variable | 0.4-1 μ M |
| <i>FS</i> TM Taq DNA Polymerase (5U/ μ l) | 0.25-0.5 μ l | 1.25-2.5U/50 μ l |
| Template DNA | variable | 10pg-1 μ g |
| Total | | 50 μ l |

1.2 Recommended PCR assay with *FS*TM PCR Buffer (Mg^{2+} free)

| Reagent | Quantity, for 50 μ l reaction | Final concentration |
|--|-----------------------------------|----------------------|
| Sterile deionized water | variable | - |
| 10X <i>FS</i> TM PCR Buffer (Mg^{2+} free) | 5 μ l | 1X |
| dNTPs (10mM each) | 1 μ l | 0.2 mM each |
| Primer I | variable | 0.4-1 μ M |
| Primer II | variable | 0.4-1 μ M |
| 25mM Mg^{2+} | variable | 1-4mM |
| <i>FS</i> TM Taq DNA Polymerase (5U/ μ l) | 0.25-0.5 μ l | 1.25-2.5U/50 μ l |
| Template DNA | variable | 10pg-1 μ g |
| Total | | 50 μ l |

Table for selection of 25 mM $MgCl_2$ solution volume in 50 μ l reaction mix:

| Final Mg^{2+} conc. | 1.0mM | 1.5mM | 2.0mM | 2.5mM | 3mM | 4mM |
|-----------------------|-----------|-----------|-----------|-----------|-----------|-----------|
| Mg^{2+} Stock | 2 μ l | 3 μ l | 4 μ l | 5 μ l | 6 μ l | 8 μ l |

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

Recommendations with Template DNA in a 50µl reaction volume

| | |
|--------------------|--------------|
| Human genomic DNA | 0.1 µg-1 µg |
| Plasmid DNA | 0.5 ng-5 ng |
| Phage DNA | 0.1 ng-10 ng |
| E.coli genomic DNA | 10 ng-100 ng |

2. Mix contents of tube. Cap tubes and centrifuge briefly to collect the contents to the bottom.

When using a thermal cycler that does not contain a heated lid, overlay the reaction mixture with 25 µl mineral oil.

3. Perform 25-35 cycles of PCR amplification as follows:

| | | |
|----------------------|---------|---------------|
| Initial Denaturation | 94°C | 3 minutes |
| 25-35 Cycles | 94°C | 30 seconds |
| | 55-68°C | 30 seconds |
| | 72°C | 10-60 seconds |
| Final Extension | 72°C | 10 minutes |

4. Incubate for an additional 2 min at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.

5. Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide staining. Use appropriate molecular weight standards.

Notes on cycling conditions

- *FS*TM Taq DNA Polymerase is for High Specificity PCR applications.
- The half-life of enzyme is >40 minutes at 95°C.

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

- The error rate of *FS*TM Taq DNA Polymerase in PCR is 2.2×10^{-5} errors per nt per cycle; the accuracy (an inverse of the error rate) an average number of correct nucleotides incorporated before making an error, is 4.5×10^{-4} (determined according to the modified method described in).
- *FS*TM Taq DNA Polymerase accepts modified nucleotides (e.g. biotin-, digoxigenin-, fluorescent-labeled nucleotides) as substrates for the DNA synthesis.
- The number of PCR cycles depends on the amount of template DNA in the reaction mix and on the expected yield of the PCR product. 25-35 cycles are usually sufficient for the majority PCR reaction. Low amounts of starting template may require 40 cycles.

Guidelines for preventing contamination of PCR reaction

During PCR more than 10 million copies of template DNA are generated. Therefore, care must be taken to avoid contamination with other templates and amplicons that may be present in the laboratory environment. General recommendations to lower the risk of contamination are as follows:

- Prepare your DNA sample, set up the PCR mixture, perform thermal cycling and analyze PCR products in separate areas.
- Set up PCR mixtures in a laminar flow cabinet equipped with an UV lamp.
- Wear fresh gloves for DNA purification and reaction set up.
- Use reagent containers dedicated for PCR. Use positive displacement pipettes, or use pipette tips with aerosol filters to prepare DNA samples and perform PCR set up.
- Always perform “no template control” (NTC) reactions to check for contamination.

Quality Control

The absence of endodeoxyribonucleases, exodeoxyribonucleases and ribonucleases is confirmed by appropriate quality tests. Functionally tested in amplification of a single-copy gene from human genomic DNA.

Endodeoxyribonuclease Assay

No detectable conversion of covalently closed circular DNA to a nicked DNA was observed after incubation of 10U *FS*TM Taq DNA Polymerase with 1 μ g pBR322 DNA for 4 hours at 37°C and 70°C.

BKL ZONE. Calle Finlandia, parcela 108, nave 5. Polígono Industrial Tecnocórdoba.
14014 Córdoba, España.

Exodeoxyribonuclease Assay

No detectable degradation of lambda DNA-HindIII fragments was observed after incubation of 10U FSTM Taq DNA Polymerase with 1 µg digested DNA for 4 hours at 37°C and 70°C.

Ribonuclease Assay

0% of the total radioactivity was released into trichloroacetic acid-soluble fraction after incubation of 10U FSTM Taq DNA Polymerase with 1 µg E.coli [3H]-RNA (40000cpm/µg) for 4 hours at 37°C and 70°C.